

HyPro technical



Hydrogen Peroxide solution designed to disinfect efficiently and very quickly all the places where hygiene is essential

Hymetec commercializes and develops products as well as technology which act effectively on any kind of infection with the aim of guaranteeing the best possible hygiene for everyone.

HyPro technical, H₂O₂ solution concentrated at 7%, provides effective and fast disinfection, spray safely

Medical Device : 0499

ISO 13485 /9001



Technical informations

Usage	Ready to use solution
Actif ingrédient	Hydrogen peroxyde 7%
GMDN Nr.	45058
Parfumes	2 flavors available: Neutral, apple
Water solubility	Total
Shelf life	17 month unopened 2 month after opening
Chemical stability	Stable from 5°C to 30°C
Relative density	1
Storage conditions	In original packaging Protect from sunlight Avoid temperature fluctuation
Packaging	Bottles of 1L (boxes of 6 bottles) Can of 5L

Efficiency

Tests were carried out on Staphylococcus aureus, Pseudomonas aeruginosa, Escherichia coli, Aspergillus Niger, Penicillium SPP, Norovirus, Bacillus stearotherophilus, ...

Applications

The HyPro technical is applicable in various fields:

hospitals, nursing homes, ambulances, dental offices, pharmaceutical, para-medical community

Warnings

- During any manipulation, wear an appropriate protection outfit (gloves, eye protection, ...)
- Check the expiry date of the product
- Read the safety sheet (available on our website: www.hymetec.net)
- Handle the product away from any source of heat

Disinfection by air is an additional tool for mechanical cleaning methods and must be in accordance with the cleaning protocol applied in the institution

Benefits

- ✓ Quick disinfection of the volumes and areas in less than 1 hour
 - ✓ Ready to use solution
 - ✓ No residue
 - ✓ No residual humidity
 - ✓ The product (H_2O_2) is ecologically inoffensive and biodegradable
- No corrosion on most of the areas.





EC Certificate Full Quality Assurance System: Certificate BE16/819942585

The management system of

Analisis s.a/n.v.

Rue de Neverlée 11
5020 Suarlée, Belgium

has been assessed and certified as meeting the requirements of

Directive 93/42/EEC
on medical devices, Annex II (excluding Section 4)

For the following products

**Chemical disinfectants intended to disinfect non-invasive
medical devices:
HyPro Medical 3% and Hypro Technical 7%**

This certificate is valid from 22/02/2017 until 28/02/2021 and
remains valid subject to satisfactory surveillance audits.
Issue 1. Certified since 22/02/2017.
Re certification audit due before 23/12/2018.

Certification is based on reports numbered BE/AMD 15/1168.QMD.

Authorised by

Pieter Weterings
Certification Manager

SGS Belgium NV, Notified Body 1639


SGS House Noorderlaan 87 2030 Antwerp Belgium
t +32 (0)3 545-48-48 f +32 (0)3 545-48-49 www.sgs.com

SGS Belgium MDD 01 0311

Page 1 of 1



Accreditation Number
005-CMS
NBN EN 17021:2011



This document is issued by the Company subject to its General Conditions of Certification Services, unless otherwise agreed, accessible at www.sgs.com/terms_and_conditions.htm. Attention is drawn to the limitations of liability, indemnification and jurisdictional issues established therein. The authenticity of this document may be verified at <http://www.sgs.com/en/Our-Company/Certified-Client-Directories/Certified-Client-Directories.aspx>. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.